

CATALOGUE #: 8CA15

PRODUCT NAME: CA15-3, standard grade

Source:	Human milk. Blood sample from the donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.
Applications:	Calibrator in immunoassays, in vitro studies.
Purification:	High-speed centrifugation, salt precipitation and delipidization.
Presentation:	Solution in 1 M NaCl, containing 0.09 % of sodium azide as preservative.
Storage:	-20°C (-15 ... -30 °C allowed)
Other information:	Avoid repeated freezing and thawing. It is recommended to aliquot the product after the first thawing.
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.