

Datasheet

Blood coagulation and Anemia • Bone Metabolism • Cardiac Markers • Fertility and Pregnancy
Gangliosides • Hormone Markers • Immunology and Serology • Infectious Diseases • Inflammation
Kidney Diseases • Metabolic Syndrome • Microbial and Plant Toxins • Miscellaneous • Neuroscience
Thyroid Diseases • Tumor Markers • Veterinary

CATALOGUE #: 8C10b

PRODUCT NAME: Calmodulin, bovine

Source: Bovine brain tissue

Application: Biochemical and immunochemical investigations.

Analysis: Purity > 95% (native PAGE in reducing conditions with EGTA).

Purification: Ion-exchange chromatography followed by gel-filtration.

Presentation: Lyophilized from protein solution with 2 mM EDTA.

It is recommended to reconstitute this product with deionized water to its initial concentration.

Storage: Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.