

CATALOGUE #: 8C10h

PRODUCT NAME: Calmodulin, human

Source: Human brain tissue. Blood samples from tissue donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Application: Biochemical and immunochemical investigations.

Analysis: Purity > 95 % (native PAGE in reducing conditions with EGTA).

Purification: Ion-exchange chromatography followed by gel-filtration.

Presentation: Lyophilized from protein solution with 2 mM EDTA.
It is recommended to reconstitute this product with deionized water to its initial concentration.

Storage: Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.