

CATALOGUE #: 8CT8

PRODUCT NAME: Canine thyroglobulin

Source:	Canine thyroid gland
Applications:	To be used as an antigen in immunoassays for determination of thyroglobulin autoantibodies in canine serum.
Analysis:	Purity > 90% (SDS-PAGE) Immunological identity confirmed by reaction with monoclonal antibody that is specific for the canine thyroglobulin.
Presentation:	Lyophilized from PBS, pH 7.4 It is recommended to reconstitute this product with deionized water to its initial concentration.
Storage:	Lyophilized and reconstituted at -20 °C (-15 ... -30 °C allowed)
Other information:	Avoid repeated freezing and thawing. It is recommended to aliquot the product after reconstitution. It is recommended to use PBS, pH 7.4, to prepare diluted samples.
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material.