

**CATALOGUE #:** 2FSH2

**PRODUCT NAME:** Monoclonal mouse anti-Follicle stimulating hormone (FSH), beta chain

<b>MAbs:</b>	F2 Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with purified natural FSH beta.
<b>Specificity:</b>	High affinity for the beta chain of FSH. No cross-reactivity with LH, TSH or HCG.
<b>MAb isotypes:</b>	<b>IgG1</b> for MAb F2
<b>Applications:</b>	EIA, Western blotting Recommended pair for sandwich immunoassay (capture-detection): F2 – F1 (Cat.# 2H8)
<b>Purification:</b>	Ion exchange chromatography
<b>Presentation:</b>	PBS, pH 7.4, 0.09 % sodium azide (NaN <sub>3</sub> )
<b>Storage:</b>	+4 °C (+2 ... +8 °C allowed)
<b>Material safety note:</b>	This product is sold <b>for research or further manufacturing use only</b> . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.