

CATALOGUE #: 3L1

PRODUCT NAME: Monoclonal mouse anti-*Listeria monocytogenes*

MAbs:	LZF7, LZH1
Specificity:	MAbs recognize <i>Listeria monocytogenes</i> outer membrane fraction and intact <i>L. monocytogenes</i> cells.
MAb isotypes:	IgG1 for MAb LZH1 IgG2a for MAb LZF7
Applications:	Indirect ELISA and Western blotting Recommended pair for sandwich immunoassay (capture-detection): LZH1 – LZF7.
Purification:	Protein G chromatography.
Presentation:	PBS, pH 7.4, 0.09 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.