

**CATALOGUE #:** PPC3

**PRODUCT NAME:** Polyclonal anti-procalcitonin (PCT)

<b>Immunogen:</b>	Human recombinant full-length procalcitonin
<b>Host animal:</b>	Goat
<b>Specificity:</b>	Specific to human PCT. Tested with human recombinant full length PCT. Cross-reactivity with other possible proteins has not been tested.
<b>Applications:</b>	Human procalcitonin immunoassay
<b>Purification:</b>	Immunoaffinity method utilizing PCT-sepharose as an affinity matrix.
<b>Presentation:</b>	PBS, pH 7.4, 0.09 % sodium azide (NaN <sub>3</sub> )
<b>Storage:</b>	+4 °C (+2 ... +8 °C allowed)
<b>Material safety note:</b>	<p>This product is sold <b>for research or further manufacturing use only</b>. Standard Laboratory Practices should be followed when handling this material.</p> <p>Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.</p>