

**CATALOGUE #:** PRP4

**PRODUCT NAME:** Polyclonal anti-canine C-reactive protein (cCRP)

**Description:** IgG fraction of goat antibodies

**Immunogen:** Canine recombinant CRP

**Host animal:** Goat

**Specificity:** Specific to canine CRP, no cross-reactivity to other serum proteins in sandwich immunoassay.

**Applications:** Canine CRP immunoassays

**Purification:** Protein G chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.