

CATALOGUE #: 8BFP

PRODUCT NAME: BNP and NT-proBNP free plasma, EDTA

Source:	<p>Pooled normal human EDTA plasma.</p> <p>Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies and HCV.</p>
Application:	<p>BNP and NT-proBNP free plasma is free of whole proBNP molecule and any proBNP fragments (BNP and NT-proBNP). It can be used for preparation of any proBNP-related standards and calibrators for immunoassays based on detection of any kinds of proBNP-related peptides.</p>
Purification:	<p>Immunoaffinity chromatography of pooled normal human plasma.</p> <p>Gel for affinity column was prepared by coupling anti-BNP and anti-NT-proBNP monoclonal antibodies with different epitope specificity with BrCN-activated Sepharose CL 4B.</p>
Presentation:	<p>Frozen liquid.</p>
Storage:	<p>-20 °C (-15 ... -30 °C allowed)</p>
Material safety note:	<p>This product is sold for research or further manufacturing use only. Standard Laboratory Practices should be followed when handling this material.</p>