

CATALOGUE #: 8FFS

PRODUCT NAME: Fatty acid binding protein (FABP) free serum

Source:	Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.
Application:	FABP free serum can be used for preparation of human FABP standards and calibrators.
Analysis:	According to ELISA testing FABP free serum does not contain more than 0.5 ng/ml of human FABP.
Purification:	Immunoaffinity chromatography. Gel for affinity column was prepared by coupling anti-FABP monoclonal antibodies (several monoclonal antibodies specific for different epitopes of FABP molecule) with BrCN-activated Sepharose CL 4B.
Presentation:	Frozen liquid.
Storage:	-20 °C (-15 ... -30 °C allowed)
Material safety note:	This product is sold for research or further manufacturing use only . Standard Laboratory Practices should be followed when handling this material.