

**CATALOGUE #:** 8MPFS

**PRODUCT NAME:** Myeloperoxidase free serum

<b>Source:</b>	Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.
<b>Application:</b>	Myeloperoxidase free serum can be used for preparation of human myeloperoxidase standards and calibrators.
<b>Analysis:</b>	No detectable amounts of human myeloperoxidase.
<b>Purification:</b>	Affinity chromatography of pooled normal human serum. Gel for affinity column was prepared by coupling anti-myeloperoxidase monoclonal antibodies (several monoclonal antibodies with different epitope specificity) with BrCN-activated Sepharose CL 4B.
<b>Presentation:</b>	Frozen liquid
<b>Storage:</b>	-20 °C (-15 ... -30 °C allowed)
<b>Material safety note:</b>	This product is sold <b>for research or further manufacturing use only</b> . Standard Laboratory Practices should be followed when handling this material.