

CATALOGUE #: 4152

PRODUCT NAME: Monoclonal mouse anti-insulin-like growth factor bindingprotein-1 (IGFBP-1 (pp12))

MAbs: G2

Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with purified IGFBP-1 from human amniotic fluid.

MAb isotypes: IgG2a for MAb G2

Applications: Indirect ELISA and Western blotting

Purification: Protein G chromatography

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.