

**CATALOGUE #:** 4IG8

**PRODUCT NAME:** Monoclonal mouse anti-insulin-like growth factor binding protein-1 (IGFBP-1)

**MAbs:** G5F8, C7B9

**Specificity:** MAbs do not cross-react with human insulin or proinsulin

**MAB isotypes:** IgG1 for MAbs C7B9, G5F8

**Applications:** Antibodies are working in ELISA and Western blotting.  
Recommended pair for sandwich ELISA (capture-detection):  
C7B9 – G5F8

**Purification:** Protein A chromatography

**Presentation:** PBS, pH 7.4, 0.04 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.