

**CATALOGUE #:** 1K9

**PRODUCT NAME:** Monoclonal mouse anti-light chains of human Ig

**MAbs:** 7A9

Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunised with mixture of light chains of human immunoglobulin.

**Specificity:** Human IgA, IgG and IgM

**MAb isotypes:** IgG2a for MAb 7A9

**Applications:** Preparation of enzyme-labeled conjugates.

**Purification:** Protein A chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.