

CATALOGUE #: 8G16-6b

PRODUCT NAME: Disialoganglioside GD1a, bovine

Source: Bovine brain tissue.

Applications: One of the main neuronal gangliosides. Useful as a differentiation marker of tumor cell growth. Can be converted to GM1 by bacterial and other sialidase from mammalian sources without detergent.

Analysis: Purity more than 98 % by HPTLC. MW 1827

Presentation: Lyophilized from water solution.

Storage: Before and after reconstitution at -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.