



# DataSheet

**CATALOGUE #:** RC4T19

**PRODUCT NAME:** Recombinant anti-cardiac troponin T (cTnT)

---

**MAbs:** RecChim406

Recombinant full-sized chimeric MAb 406 expressed in mammalian cells. Original wild type variable domains of MAb 406 (Cat.# 4T19) and human IgG1 constant domains are combined in this construct.

**Specificity:** Human cardiac troponin T

**Applications:** Human cardiac troponin T sandwich immunoassays

**Purification:** Protein A chromatography

**Presentation:** PBS, pH 7.4, 0.09 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research or further manufacturing use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligent 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)